

STP 310

Revolutionary Infectious Substance Packaging

The STP 310 meets the requirements outlined in IATA Packing Instruction 602 for safely shipping infectious substances. The STP 310 also meets the requirements for IATA Packing Instruction 650 for shipping large volumes of diagnostic specimens (up to 4 Liters)

Holds more

243 - 2 mL tubes
2 x 500 mL Blood bags
32 x 10 mL Evacuated Blood Collection tubes
many other configurations approved call for listing

Costs less


Costs from 50% to 90% Less per specimen than conventional Packaging

Accepting Orders for December Delivery

Disclaimer:

This package meets the test criteria outlined in federal regulations for the safe transport of infectious substances if used as directed. Shipper must use the package system as directed to ensure compliance.

Manufacturer/Agent/Distributor will not be liable if the package system is altered or not used according to these instructions. Federal regulations require that before shipping dangerous goods, the Shipper must be conversant and comply with applicable shipping regulations.



10807-184 Street
Edmonton Alberta
Canada,
T5S-5T1

Tel: (780) 486-0211
1 800 814-7484
Fax: (780) 486-0235
<http://www.saftpak.com>