

STP-350

Insulated Infectious Substance Packaging

Certified EPS Insulated Packaging for Infectious Substances.

Meets I.C.A.O Packing Instruction #602 (Infectious Substances)

Meets I.C.A.O Packing Instruction #904 (Dry Ice)

Meets I.A.T.A Packing Instruction #650 (Diagnostic Specimens)

Meets Transport Canada TC-125/1A Requirements



Suitable for up to 4 Litres on Passenger and Cargo Aircraft using I.A.T.A. special provision A81

Internal secondary containers available separately

Temperature remained between 0 and -44 degrees celcius for 70 hours

TC-125-1A/8-14

UN 4H CLASS 6.2/97
CAN/8-14 SAF-T-PAK



SAFTPAK[®]
INC.

Edmonton Alberta
Canada,

Tel: (403) 486-0211
1 800 814-7484
Fax: (403) 486-0235
<http://www.saftpak.com>